


Characterization of the ongoing electrophysiological activity of the planarian *Schmidtea mediterranea*



Julian Keil, Lukas Lang, Maren Eberle, Jannes Freiberg, Christian Kaernbach
Christian-Albrechts-Universität

www.biopsych.uni-kiel.de | keil@psychologie.uni-kiel.de | @drjuliankeil

https://github.com/juliankeil/Slides/blob/master/Talk_Teap_2022.pdf

Keep it simple!


What is the simplest possible model system?

- Bilateral organization
- Visual processing
- Cognitive processes

Schmidtea mediterranea as an animal model, to study functional neural mechanisms in simple nervous systems

Close to the common ancestor of bees and humans

Fraguas et al., 2012




Previous work

Christian-Albrechts-Universität zu Kiel

Invasive recordings from the cephalic ganglia


- Animal cooled to 5-7° C
- Electrode inserted 0.2 – 0.3 mm
- Ongoing recordings for 20-30 minutes until the death of the animals



Ongoing oscillations?

Christian-Albrechts-Universität zu Kiel

- 2 animals
 - Warmed to $\sim 10^\circ \text{C}$
- No separation of ongoing activity and muscle potentials
- FFT analysis
 - 0.1 – 10 Hz
 - $1/f$ -characteristic above 0.5 Hz



Optimize recordings


Christian-Albrechts-Universität zu Kiel

Goals:

- Less cooling
- No harming
- Simple setup


Approach

- Human EEG amplifier
- Electroretinography wire electrodes
- Fixate the animal in agarose gel
- State-dependent recordings



Optimize recordings - Setup


- Pipette animal on specimen slide
 - Slide is cooled on crushed ice
- Pipette agarose gel on animal
 - Gel max. 50° C warm
- Slide hardened gel with animal onto wire electrodes
- Ongoing recordings for 10 minutes
 - Darkness = 0 lux
 - Light = 40000 lux
- Slide animals into petri dish for observation and recovery



First results

Christian-Albrechts-Universität zu Kiel

- N = 17 darkness, N = 20 light
- 10 minutes ongoing recording
 - 0.1 – 20 Hz FIR filter
 - 10s segments, exclude muscle activity
 - 0.5 – 20 Hz FFT
- Between-animals independent-samples t-test with cluster correction for multiple comparisons
 - 1/f-characteristic during darkness
 - Broadband increase in power between 3 and 20 Hz during light stimulation



Mean light

Individual light

Mean darkness

Individual darkness

What's next?

Current state:

- Successful recording of ongoing neural activity without harming or killing the animal
- First evidence of functional changes of ongoing neural activity due to light stimulation

Next step:

- Optimize recording environment to exclude noise
- Introduce event-based experiments



Thank you!



Christian-Albrechts-Universität zu Kiel

Lukas Lang, Maren Eberle, Jannes Freiberg,
Christian Kaernbach




References

- van den Heuvel, M. P., Bullmore, E. T., & Sporns, O. (2016). Comparative Connectomics. *Trends in Cognitive Sciences*, 20(5), 345–361. <http://doi.org/10.1016/j.tics.2016.03.001>
- Popov, T., & Szyszka, P. (2020). Alpha oscillations govern interhemispheric spike timing coordination in the honey bee brain. *Proceedings of the Royal Society B: Biological Sciences*, 287(1921), 20200115. <http://doi.org/10.1098/rspb.2020.0115>
- Fraguas, S., Barberán, S., Ibarra, B., Stöger, L., & Cebri, F. (2012). Regeneration of neuronal cell types in *Schmidtea mediterranea*: an immunohistochemical and expression study. *The International Journal of Developmental Biology*, 56(1-2-3), 143–153. <http://doi.org/10.1387/ijdb.113428sf>
- Aoki, R., Wake, H., Sasaki, H., & Agata, K. (2009). Recording and spectrum analysis of the planarian electroencephalogram. *Neuroscience*, 159(2), 908–914. <http://doi.org/10.1016/j.neuroscience.2008.11.011>


Optimize recordings -Setup



Human EEG amplifier



State-dependent recordings



Electroretinography wire electrodes